# 甲基化芯片数据分析

### 主要内容

- ➤ 从GEO下载数据
- ▶ 质控
- > 归一化
- > 差异甲基化位点
- > 差异甲基化区域
- > 注释

#### 甲基化芯片

Illumina Human Methylation 450K BeadChip芯片可 检测人全基因组近450,000个甲基化位点,具有单碱 基的分辨率。全面覆盖了96%的CpG岛,并根据需 求加入了CpG岛以外的CpG位点、人类干细胞非 CpG甲基化位点、正常组织与肿瘤(多种癌症)组 织差异甲基化位点、编码区以外的CpG岛、miRNA 启动子区域和GWAS疾病相关区域的位点,同时覆 盖了Human Methylation27 BeadChip的90%的位 点。

#### 芯片的值

- A探针(非甲基化)的数目U
- B探针(甲基化)的数目M
- β值或者m值

$$\beta = \frac{M}{(M+U+100)}$$
 m=log  $(\frac{\beta}{1-\beta})$ 

- β值反映了能够和给定被甲基化的序列匹配的寡核苷酸的 比率,序列中的甲基化率
- M值可以消除探针不同而造成的影响

### 数据下载


1. **GEO** 

http://www.ncbi.nlm.nih.gov/geo/

2. TCGA、EBI等其它网站


### 质控

#### Beta





### 质控

MDS 1000 most variable positions



### 归一化






# 差异甲基化位点

| Id            | intercept    | f           | pval     | qval        |
|---------------|--------------|-------------|----------|-------------|
| cg19229544    | -5.819817249 | 9137.68304  | 7.18E-08 | 0.013840388 |
| cg13706769    | 3.192061949  | 7427.708505 | 1.09E-07 | 0.013840388 |
| ch.1.1262625R | -16.60962605 | 6836.276705 | 1.28E-07 | 0.013840388 |
| cg01972335    | 2.144506398  | 5554.829221 | 1.94E-07 | 0.01571845  |
| cg19506433    | -0.65297139  | 4201.708582 | 3.39E-07 | 0.02196957  |
| cg05945291    | -16.60962605 | 2652.551359 | 8.51E-07 | 0.033324073 |
| cg21081097    | -6.283964656 | 2621.428714 | 8.71E-07 | 0.033324073 |
| cg08094135    | 2.349302795  | 2547.055433 | 9.22E-07 | 0.033324073 |
| cg25010146    | -5.218557921 | 2449.915749 | 9.97E-07 | 0.033324073 |

# 热图



# 差异甲基化区域

| chr   | start         | end           | value                | area            | cluster | indexSt<br>art | indexE<br>nd | L  | cluster<br>L | p.value | fwer |
|-------|---------------|---------------|----------------------|-----------------|---------|----------------|--------------|----|--------------|---------|------|
| chr1  | 2525719<br>1  | 2525823<br>6  | 0.43417<br>0033      | 8.68340<br>0661 | 4866    | 11732          | 11751        | 20 | 33           | 0       | 0    |
| chr11 | 2397201       | 2397977       | 0.42381<br>6783      | 8.05251<br>8884 | 31573   | 274030         | 274048       | 19 | 37           | 0       | 0    |
| chr6  | 3193837<br>2  | 3193911<br>2  | -<br>0.43463<br>8412 | 7.38885<br>3009 | 163981  | 163875         | 163891       | 17 | 109          | 0       | 0    |
| chr6  | 1498059<br>95 | 1498067<br>32 | -<br>0.55850<br>1181 | 6.14351<br>2994 | 170492  | 181364         | 181374       | 11 | 11           | 0       | 0    |
| chr2  | 2092242<br>25 | 2092251<br>83 | -<br>0.61795<br>2848 | 3.70771<br>7086 | 115900  | 73375          | 73380        | 6  | 8            | 0       | 0    |

#### 差异甲基化区域

#### value--average smooth methylation difference within block

area--block area estimate (abs(value) \* length)

cluster--id of cluster blockgroup within which block occurs

indexStart--index of first cluster in block

indexEnd--index of last cluster in block

L--number of clusters in block

clusterL--number of probes in block

p.value--permutation p.value based on difference conditioned on length

fwer--family-wise error rate estimate based on difference conditioned on length

p.valueArea--permutation p.value based on area

fwerArea--family-wise error rate estimate based on area

## 差异甲基化区域注释

| Chr   | Start     | End       | Ref | Alt | Func           | Gene     |
|-------|-----------|-----------|-----|-----|----------------|----------|
| chr1  | 25257191  | 25258236  | 0   | 0   | intronic       | RUNX3    |
| chr11 | 2397201   | 2397977   | 0   | 0   | ncRNA_intronic | CD81-AS1 |
| chr6  | 31938372  | 31939112  | 0   | 0   | exonic         | DXO      |
| chr6  | 149805995 | 149806732 | 0   | 0   | UTR5           | ZC3H12D  |
| chr2  | 209224225 | 209225183 | 0   | 0   | downstream     | PIKFYVE  |
| chrY  | 6778641   | 6778641   | 0   | 0   | upstream       | TBL1Y    |

注释网站: http://wannovar.usc.edu/